



4621 East 75th Terr. Kansas City, MO 64132

Phone (816) 300 2828 | Fax (816) 300 0864

e-mail [livers@liversbronze.com](mailto:livers@liversbronze.com)

## 1. PRODUCT NAME

Handrailing system bearing the name BELMONT

## 2. MANUFACTURER

Livers Bronze exclusively manufactures the handrailing system as described on this sheet.

## 3. BASIC USE

Livers Bronze BELMONT railing system is a pre-engineered component based railing system, shipped to the jobsite ready to install.

**LIMITATIONS:** BELMONT railing system is only limited by the materials from which it is fabricated.

Pre-engineered railing should be assembled using the components listed:

42" high guardrail post for tempered glass infill applications. Shall be powder coated cold rolled steel or type 304 #4 stainless steel (Belmont-SS) spaced at 4' - 0" O.C. maximum.

42" high guardrail post for tempered glass infill applications w/side mounted handrail at 34" above stair nosing. Shall be powder coated cold rolled steel or type 304 #4 stainless steel (Belmont-SS) spaced at 4' - 0" O.C. maximum.

**TOP RAIL:** 1.5' round steel, stainless steel, aluminum or bronze. (2" round minimum for wood top rail.)

**ASSIST RAIL:** 1.5" round steel, stainless steel, aluminum or bronze.

**FINISHES:** A wide selection of colors and textures in powder coat is available for post finishes. Colors are selected from manufacturers standard color charts. Standard circular satin and polished finishes are available for stainless steel and bronze railings. Anodized finishes are not available.

**GLASS:** 3/8" tempered glass infill furnished by Livers Bronze for field attaching to BELMONT post.

## 4. TECHNICAL DATA

The BELMONT railing system is designed to meet performance criteria as outlined in IBC 2012.

## 5. INSTALLATION

Livers Bronze BELMONT railing system must be attached with fasteners and to surfaces that are approved or selected by the manufacturer, and must be detailed on all shop drawings.

## 6. AVAILABILITY AND COST

**AVAILABILITY:** BELMONT railing system is subject to current delivery schedules. Each project must be reviewed before scheduling is available.

**COST:** Varies according to quantity and material. Further information may be obtained from Livers Bronze or one of its authorized representatives.

## 7. GUARANTEE

Livers Bronze warrants that its BELMONT railing system is of sound quality, conforming to standards as stated in part 4 hereof and is reasonably free of defects in material and workmanship. Inasmuch as Livers Bronze has no control over the installation of its products, no other warranty is expressed or implied. Livers Bronze does not warrant against glass breakage.

## 8. MAINTENANCE

The surface of Livers Bronze material can be kept clean by the use of ordinary soaps and mild cleansers. Abrasive cleaners should not be used.

## 9. TECHNICAL SERVICES

Technical services are rendered direct from Livers Bronze in Kansas City, Missouri.